ELECTRONIC KIDS WHEELCHAIR



FEATURES

Kids wheelchair for every terrain

The Forest Kids is an electronic wheelchair for children, used for indoor and outdoor and ideal in overcoming obstacles. The 4 wheels are always in contact with the ground even on rough surfaces up to 120 mm.



orange



inox













Electronic tilting of seat and back

Height adjustable pushing handles

Armrests adjustable in height







Legrests adjustable in width





Stability axle





















Wide range of accessories and options.

2x12V-75Ah 2x220W 2x350W

10°/17%





*iPortal is compatible with all the electronic wheelchairs of Vermeiren equiped with a DX2 Basic, DX2 LCD or special controller.

























600

15°

200-250

530

330-390

350

-5°/55°

50