

# EZ1

COMFORTABLE  
RELAXING CHAIR



# Technical information EZ1

total width	seatwidth	height armrest	seat height	total height	total length	seat depth	back height	adjustment seat	adjustment back	adjustment legrest	weight user
60	39										
63	42										
66	45										
69	48	19-29	49-62	102	118,5	43-52	50-60	-5°/+15°	30°	19°/113°	max.150
72	51										
75	54										
78	57										

Deviations in sizes and weight (1,5 cm and 1,5 kg) are possible


## Possibilities




Adjustable headsupport in all directions



Removable & adjustable armrests / seat



Fully adjustable - collapsible and removable feet / legrests



## General options



Side support and table



Type 30



## N.V. Vermeiren N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

Tel.: +32 (0) 3 620 20 20  
Fax: +32 (0) 3 666 48 94

[www.vermeiren.com](http://www.vermeiren.com) - [info@vermeiren.be](mailto:info@vermeiren.be)



# VERMEIREN

Opening up new horizons in mobility

For more information, contact your local dealer



Subject to technical alterations