

Bobby

Transfer chair



Transfer wheelchairs need to be light, easy to handle and foldable in mere seconds.

With its removable footrests, foldable backrest and two pull straps to carry your Bobby with you during transport, it fits the description of a transfer wheelchair perfectly. On top of that, it comes with a safety belt, brakes and adjustable footrests.



VERMEIREN
we care for you

Comfact foldable transfer chair

The **Bobby** is your ideal compact, foldable transfer chair. Light and easy to fold in seconds.

- Foldable
- Two manual brakes
- Padded armrests
- Angle-adjustable footplates
- Adjustable heelstraps

Either on holiday, in a hospital or airport for a transfer from point A to point B, the Bobby chair is easy to fold and to store and requires little maintenance. Also available with 24" rear wheels.



Foldable back with pouch



Tipping piece at the back



Optional removable side cover



Foldable for easy transport

Your guarantee of quality

Regarding quality, the **Bobby** complies with the highest norms, laid down by certified test institutes in most European countries.

Technical data

■ Total width	575 / 635 mm
■ Total height	940 mm
■ Total length	950 mm
■ Folded width	270 mm
■ Seat width	420 / 480 mm
■ Height armrests	220 mm

■ Seat height	500 mm
■ Height backrest	380 mm
■ Seat depth	420 mm
■ Weight	> 9,4 kg
■ Maximum user weight	115 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you