



376

Overbed Table

The bed table 376 is height adjustable, can be tilted and has casters for easy mobility. Ideal for home or hospital use, in combination with a bed, sofa or a recliner. The frame is finished with long-lasting powder coated paint and the table itself can be tilted and locked in both directions. Ideal for reading while the smaller, stationary part is perfect to place drinks, food or electronic devices.

■ Total length	920 mm
■ Total width	470 mm
■ Adjustable height	610 - 950 mm
■ Underframe height	120 mm
■ Tabletop	765 x 410 mm
■ Total weight	9,75 kg
■ Max. belasting	15 kg
■ Wheels	ø 50 mm
■ Colour(s)	frame : grey - table: oxide frame: beige - table: beech
■ Additional information	CE

Measurement tolerances: ± 15 mm / 1,5kg. We reserve the right to introduce technical changes. Models shown only for illustration purposes. Certain features may be optional.

