

Inovys II



- Adjustment of the armrests in length and height
- Adjustment of the seating depth
- Adjustment of the footrests in width
- Possible to lower the chair without changing its original wheels



Tipping pieces



L58 - Adjustable headrest



Headrest adjustable in angle, height and depth with rotation and translation

Armrests can be removed. Adjustable in height, length and width.

Pushing handles adjustable in height



Memory foam seat cushion




Optional L04 Side supports



Optional B12 - Table

Wide range of accessories and options.



620	390-440	190-280	420-480	1030	500-560	0°-20°	35°	19°-113°	420-480	1100	135
670	400-490	190-280	420-480	1030	500-560	0°-20°	35°	19°-113°	420-480	1100	135
720	390-440	190-280	420-480	1030	500-560	0°-20°	35°	19°-113°	420-480	1100	135

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.