

1030 - 1035

Electrical and manual adjustable bed



1030

1035

N.v. Vermeiren N.v.

Vermeirenplein 1-15

B-2920 Kalmthout

Belgium

Tel.: 00 32 3-620 20 20

Fax: 00 32 3-666 48 94

info@vermeiren.be

www.vermeiren.com



1030

Electrical high/low bed

electrical head- and footrest




| | |
|-------------|---------|
| 1030 | |
| headrest | 16,13kg |
| footrest | 15,34kg |
| total | 31,47kg |
| slats | 4,54kg |
| siderails | 6,91kg |
| selfraiser | 6,3kg |
| front panel | 16,13kg |
| rear panel | 16,13kg |



Clear markings allow easy and correct plug-in.


The self-raising device was placed outside the useful area and creates as such more available space for the mattress. Mattresses up to 2 meters in length can be used.



1035

Standard bed

hydraulic headrest and manual footrest adjustment



1,80m x 0,52m.
Compact and easy to fold.

| | |
|-------------|---------|
| 1035 | |
| headrest | 15,87kg |
| footrest | 11,56kg |
| Total | 27,43kg |
| slats | 4,54kg |
| siderails | 6,91kg |
| selfraiser | 6,30kg |
| front panel | 10,8kg |
| rear panel | 10kg |

Maximum weight: 170kg
Subject to technical alterations

