

Technical information Inovys II

total width	seat width	folded width	seat inclination	back inclination	height armrests	height seat	total height	total length	seat depth
630	400	-	0°-21°	0°-54°	220-320	460-530	1020-1200	1010-1050	460-520
680	450	-	0°-21°	0°-54°	220-320	460-530	1020-1200	1010-1050	460-520
730	500	-	0°-21°	0°-54°	220-320	460-530	1020-1200	1010-1050	460-520

Complies to crash test ISO 7176-19

Deviations in sizes and weight (15 mm and 1,5 kg) are possible

height back	maximum inclination	front wheels	rear wheels	max. weight	users weight	TUV ISO 9001	CE	Crash Test Approved
550-585	15°	200	16" - 22" - 24"	36	135			
550-585	15°	200	16" - 22" - 24"	36	135			
550-585	15°	200	16" - 22" - 24"	36	135			

General options

- Table
- Safety belt
- Spreader
- Calf support with dartex padding
- Cannister holder

Colors

C55
inox

* Other colors upon request.

The wheelchair must be adapted to you, not vice versa

Our basic principle is: the wheelchair must be adapted to you, not vice versa. This is why our wheelchairs have a modular design. All types of wheelchairs can be adapted, have numerous accessories and all options are mutually interchangeable. This is the result of our continuous research and more than 55 years of experience.

Your guarantee of quality

Regarding quality, the Inovys II complies with the highest norms laid down by certified test institutes in most European countries..

- Adjustment of the armrests in length and height
- Adjustment of seat depth
- Adjustment of the legrests in width
- Possible to lower the height of your chair without changing to smaller wheels

Inovys II & Inovys II max. comfort



N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be



VERMEIREN
Opening up new horizons in mobility

For more information, contact your local dealer



Information for the visually impaired is available upon request. - R.E.: N.V. Vermeiren N.V., Vermeirenplein 1/15 - 2920 Kalmthout - Belgium - 12/2011



model shown is for illustration only



Opening up new horizons in mobility



Push handles adjustable in height



Metal push bar in one piece



Removable armrests. Adjustable in height and depth. Adjustable seat depth.



New wooden removable table



Easy to disassemble

Inovys II

Wheelchair equipped with 24" Krypton rear wheels which can be removed with a Quickrelease system. Front wheels are 200 mm x 35 mm. Armrests can be adjusted in height and width. Headrest can be adjusted in height, depth and width.



The new comfortable Inovys II brings an optimum comfort to the user. It's equipped with a back and a seat in Dartex for a perfect seating comfort. A new attachment system for the backrest allowing adjustments in depth, height and angle. The seat is maintained by elastic straps for more flexibility and comfort making it possible to avoid pressure points.

The Inovys II can be totally inclined, is easy to fold and even easier to disassemble for transport.

NEW

very solid attachment of the legrests for a long life usage in residence or community.

Inovys II maxi confort



Comfortable headrest adjustable to your personal needs

Standard equipped with :

- Guttered armrests
- Comfort headrest L58C
- Comfort leg rests BZ8C (dartex paddings)



Hemiplegic directional armrest - B66O



L04 - Side supports



T30 version with 16" rear wheels



Feet rubbers for easy tipping