

# Carpo 3 & 4

## FEATURES

Speed and longlasting pleasure



Electronically programmable

brown



Adjustable swivelseat with removable headrest and adjustable armrests

### Model Carpo 3 Limited edition

Powerful outdoor comfortable 3-wheel scooter. Electronic programmable and a reliable partner for outdoor trips.

optional  
**80A batteries**  
available!



Lighting front and rear and indicators

Bumper

Inflatable tires

Extra strong shoppingbasket

### Model Carpo 4 Limited edition

optional  
**80A batteries**  
available!

Powerful and comfortable 4 wheel outdoor scooter. The Carpo 4 Limited edition is electronically programmable and standard equipped with a bright LCD with various indications.



Lighting front and rear and indicators

Bumper

Inflatable tires

Adjustable seat

Also available as  
**Model Carpo 4**  
in champagne color



champagne

available  
**XD version**  
For hilly regions



km/h	km			
C3 15	45	9%/15%	100	3030
C4 15	45	9%/15%	100	3160

mm	mm	mm	mm	kg	kg	kg	Ah	W
C3 660	1200	1400	4.00-4/3.00-8	127	150	2 x 12V / 75 Ah	750 W	
C4 660	1350	1450	3.00-8	131	150	2 x 12V / 75 Ah	750 W	

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

CARPO 3 - 4

