

Quava

Rollator

The Quava is a modern, lightweight, foldable rollator with multifunctional brakes. Ideal for taking along in the car, during a trip or when shopping. Due to the wide, padded back belt you can, with the brakes locked, sit down and rest.

■ Total width	655 mm
■ Total length	710 mm
■ Seat height	540 mm
■ Seat width	450 mm
■ Push handles	830 - 955 mm
■ Wheels (front/rear)	200 x 35 mm
■ Multifunctional brakes with lock	
■ Weight	6,2 kg
■ Max. users weight	136 kg
■ Max. load basket	5 kg
■ Colour(s)	Carbon grey or red
■ Additional information	CE

Measurement tolerances: ± 15 mm / 1,5 kg. We reserve the right to introduce technical changes. Models shown only for illustration purposes. Certain features may be optional.

