

Twin-light

R O L L A T O R



total width	570 mm
seat depth	200 mm
seat width	390 mm
seat height	540 mm
total height	760 - 950 mm
total length	570 mm
weight	3 kg
max. users weight	130 kg
wheels	113 mm
colour	carbon grey

Measurement tolerances: ± 1.5 mm / 1,5kg.
We reserve the right to introduce technical changes.



Model shown only for illustration purposes. Certain features may be optional.

N.V. Vermeiren N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium

Tel.: +32 (0)3 620 20 20 • Fax: +32 (0)3 666 48 94 • info@vermeiren.be • website: www.vermeiren.com



Opening up new horizons in mobility