


Material			
No.	Plan	Article	Quantity
1	XP02835 (High 100)	2994644	1
	XP02836 (High 120)	2994645	1
	XP03189 (Extension High 100)	2996888	1
	XP03190 (Extension High 120)	2996889	1
2	SANMED C-460/5202 Blue	2994651	4,1m
	SANMED C-460/5202 Blue (Extension)	2994651	4,5m
3*	Slide fastener	1995604	4,1m
	Slide fastener (Extension)	1995604	4,5m
4*	Slide fastener head	1995605	2
5*	Sewing thread	W200016	20m
	Sewing thread (Extension)	W200016	21m

Final dimensions of mattress: 1950x870x100mm or 1950x870x120mm  
 Final dimensions of mattress with extension: 2150x870x100mm or 2150x870x120mm /  
 Wymiary końcowe materaca: 1950x870x100mm lub 1950x870x120mm  
 Wymiary końcowe materaca z przedłużką: 2150x870x100mm lub 2150x870x120mm

EP00468 / A

\* not depicted / nie ujęto na rysunku

MATTRESS FOR LUNA BEDS				Plan No.
				EP00468/ A
Scale: 1:15	Drafted	Verified	Confirmed	
Date	28/10/2019	04/02/2021	04/02/2021	
Initials	AZ	EB	TK	