



Use glue
WEICON AN 302-43
(H900009)



2994874

SE04(E) CHIN CONTROL TYPE LINX				Plan No.
Sheet 1/2				E001143 / F
Scale: 3:10	Drafted	Verified	Confirmed	
Date	11/06/2020	13/05/2022	13/05/2022	
Initials	SD	SD	PJ	

E001143 / F


Material			
No.	Plan	Article	Quantity
1	X005419	1910801	3
2	X005951	2078399	4
3	Y023795	2075399 (R)	1
	Y023796	2075499 (L)	1
4	Y023798	2075599 (R)	1
	Y023799	2075699 (L)	1
5	Y021921	2020154	1
6	Y021922	2020099	1
7	Y023801	2075799	1
8*	Y023959	2078499	1
9*	Y023960	2078599	1
10	M4x10 DIN84	1001464	2
11	M5x12 DIN912	1001489	2
12	M8x40 DIN912	1001500	1
13	D40 M8x20	1001580	1
14	M8 DIN985	1001682	3
15	M6x45 DIN912	1003109	1
16	D30 M6x10	1003658	1
17	M5x10 ISO7380	1007063	2
18	M6 DIN986	1008492	1
19	DLX CR400-A	1913999	1
20A*	M8x20 ISO7380	1992349	2
20B**	M8x25 ISO7380	1990501	2
21	PB30 Blue	1997690	2
22	PB30 Green	1997691	1
23**	TA38	1916865	1
24*	M5x6 DIN912	1997822	1
25***	Linx 2m	1912412	2

Material			
No.	Plan	Article	Quantity
26	X004049	1902940	1
27**	M8x10 DIN916	1008942	1
28**	Y024073	2081554	1
29**	M6x10 ISO7380	1900641	4
30**	Y024076	2081600	1
31**	X002762	1554454	1
32**	Y005639	1555100	1
33*	M8 DIN125	2994415	2
34	XC00874	1917371	1
35**	X003826	1009273	3
36	X005984	2082600	1

* Manual version SE04
 ** Electric version SE04E
 *** Not depicted
 (34) = option SE63

2994874

E001143 / F

SE04(E) CHIN CONTROL TYPE LINX				Plan No.
Sheet 2/2				E001143 / F
Scale: 3:10	Drafted	Verified	Confirmed	
Date	11/06/2020	13/05/2022	13/05/2022	
Initials	SD	SD	PJ	